



Mageshape machine+Chair 2 in 1

MAGESHAPE MIX



Muscle re-build of incontinent patients
Body shaping and muscle building



TIRED OF THE CONSTANT URGE TO GO?

Pelvic floor muscles are a layer of muscles that support the pelvic organs. Strong pelvic floor muscles provide control over the bladder and bowel. When these muscles are weakened, they do not fully support the internal organs, which may lead to difficulty controlling the release of urine.

ARE YOU ?

Suffering from these afflictions
Growing Need for Intimacy?



Dropped Wellbeing

Dissonance incouple's life



Coughing or Sneezing

Leaking urinewhile laughing




Decreased PelvicFloor
Function


Early AgingMenopause

Traditional exercise

KEGEL EXERCISE AND PELVIC EXERCISE



1. Preventing urinary incontinence and vaginal prolapse
2. Preventing uterine prolapse
3. Strengthening the pubococcygeal muscle



1. Male Prostate Pain
2. Benign prostatic hyperplasia
3. Prostatitis Magna Erectile Dysfunction

New tech for improving urinary incontinence

-Mageshape MIX

What's Mageshape MIX?

A breakthrough treatment for incontinence and confidence. This unique technology revolutionizes women's intimate health and wellness by providing those suffering from incontinence with a completely non-invasive option.



PRODUCT FEATURES



Improve sexual sensitivity
Sensitivity and function



Prevention and improvement
Stress urinary incontinence




Postpartum recovery
Improve vaginal health



Strengthen the pubic tail
Skeletal muscle (PC muscle)



Increase blood in the pelvic floor
Flow and nerve function



Increase the vagina after menopause
Wall thickness and lubricity

- ◆ Improve urine leakage and frequent urination.
- ◆ One course of treatment = 20,000 Kegel exercises
- ◆ Health care of the pelvic cavity
- ◆ Sex life is no longer affected
- ◆ No need to take off pants, expensive invasive treatment
- ◆ Both men and women can use



How Does it Work?


Utilize high intensity focused electromagnetic techniques to stimulate the deep muscles of the pelvic floor and restore neuromuscular control. A single focused electromagnetic wave energy brings the greatest contraction on thousands of pelvic floor muscles, which is a very important muscle for incontinent patients. It makes the pelvic floor muscles perform hyper-contraction exercises (20,000 pelvic floor muscle contractions every 30 minutes).




HIFEM stimulates deep layers with third-generation happiness ch




90% of people cannot exercise deep muscles properly



Unable to know whether the pelvic floor muscles are contracting correctly



The contraction frequency, intensity and amplitude are limited



100% contraction of deep pelvic floor muscles




Magnetic energy accurately focuses the pelvic floor muscle penetration depth upto 10 cm



Three-stage treatment mode is effective for training and strengthening



The strength and frequency of pelvic floor muscle contraction is limited



Ultra-limit contraction of pelvic floor muscles and maintain for a few seconds

Why choose Mageshape Mix

➤ **Non-invasive**

Mageshape MIX treatment is non-invasive and does not require needles at all. HIFEM Energy stimulates the entire pelvic floor area and helps with restoring neuromuscular control.

➤ **Improve quality of life**


95% of patients have reported significant improvement in the quality of life.

➤ **Comfortable**

Mageshape MIX treatment is comfortable. You will experience a little tingling sensation throughout the 30 minutes at your pelvic area. No numbing cream is required.

➤ **Effective**

100% contraction of deep pelvic floor muscles, magnetic energy accurately focuses on the pelvic floor muscles with a penetration depth of 10 cm



30 = 20.000

Mins Treatment

Kegel exercise

30 MINS TREATMENT
= 20.000 KEGEL EXERCISE





**2 treatment handles for
targeted booy parts training**



Belly




Buttocks



Thighs



Belly




Legs



Arms

HOW IT WORKS?



- » It use HIEMT (High-Intensity Focused Electromagnetic Wave) technology can target motor nerve cells to stimulate muscle contraction with high intensity, which deeply reshape the internal structure of muscles, strengthen the muscles, increase muscle density and volume.
- » Muscle extreme contraction requires a lot of energy supply, then the fat cells next to the muscle will be consumed, so it can help reduce the fat effectively.
- » It makes the body lines more prominent, "mermaid line" and "The vest line" clearer.
- » It can be said that "you can do exercise while lying down." In 30mins treatment, it equals 20000 workouts.

MINS TREATMENT

WORKOUTS

30 = 20.000

16% ↑ Average increase in muscle mass

19% ↓ Average fat reduction



Operation System



Sport Mode

S1 S2

TIME FLUENCE

45 09 45 09

Min % Min %

MAIN < 06 > Hz

TIME < 04 > S

SUB < 05 > Hz

TIME < 08 > S

Gentle Mode

S1 S2

TIME FLUENCE

45 09 45 09

Min % Min %

MAIN < 06 > Hz

TIME < 04 > S

SUB < 05 > Hz

TIME < 08 > S

Professional Mode

S1 S2

TIME FLUENCE

45 09 45 09

Min % Min %

MAIN < 06 > Hz

TIME < 04 > S

SUB < 05 > Hz

TIME < 08 > S

Operation System



Sport Mode

TIME 45 Min

FLUENCE 06 %

Power button, minus, plus, and activity icons (running, walking, cycling) are visible.

Gentle Mode

TIME 45 Min

FLUENCE 06 %

Power button, minus, plus, and a menu icon are visible.

Professional Mode

TIME 45 Min

FLUENCE 06 %

MAIN < 06 > Hz, TIME < 04 > S, SUB < 05 > Hz, TIME < 08 > S

Power button, minus, plus, and sub-timer settings are visible.

Operation System



Sport Mode

S1 S3

TIME FLUENCE

45 09 45 09

Min % Min %

MAIN < 06 > Hz

TIME < 04 > S

SUB < 05 > Hz

TIME < 08 > S

Gentle Mode

S1 S3

TIME FLUENCE

45 09 45 09

Min % Min %

MAIN < 06 > Hz

TIME < 04 > S

SUB < 05 > Hz

TIME < 08 > S

Professional Mode

S1 S3

TIME FLUENCE

45 09 45 09

Min % Min %


MAIN < 06 > Hz

TIME < 04 > S

SUB < 05 > Hz

TIME < 08 > S

Android system



Beauticians administration • Patients administration

Clients information storage

Treatment parameter date memory system

Rental model

Specification

Magnetic wave(energy)	0-7tesla
Voltage	220V/110V
Out put power	3200W
Frequency	F1:1-10HZ F2:1-100HZ
Pulse width	300us
Mode	Model-I(Sport mode) Model-II(Gentle mode) model-III(Professional mode)
Screen	10.4 inch
Machine size	830*715*780mm

Frequently Asked Questions

I have a low tolerance for pain! Is it painful? **A**

There's no pain at all. The experience is a tingling sensation with your pelvic floor muscles' contracting during the procedure. You may resume your daily activities immediately after the treatment.

How long does the treatment take? **B**

A typical treatment would take 30-minutes. During the treatment, you will feel a slight tingling with no pain at all.

How many sessions do I need? **C**

Always consult your provider for a tailored treatment plan. For optimal results, it's recommended to have a minimum 6 sessions, scheduled twice a week.

How long will the result last? **D**

Our patients reported the results last up to 6 months. With regular kegel exercises, it can help to sustain the results past the 6-months mark!

How about downtime and side effects? **E**

There is absolutely no downtime and side effects from the treatment. After the session, you can immediately return to your daily activities.